

Improve Better Patient Outcomes with Longer Lasting Restorations.

Increase office productivity and save time and material costs



OFFER VALID: OCTOBER 1, 2023 - DECEMBER 31, 2023

Specialty Carbides, Diamonds, Polishers and All Procedural Kits

BUY 4 GET 1 AT NO CHARGE**

CARBIDES: GREAT WHITE® ULTRA, ENDOGUIDE® AND ENDO SAFE END
DIAMONDS: PIRANHA®, REVELATION®, TDA AND GREAT WHITE®Z DIAMONDS
POLISHERS: JAZZ® POLISHERS

CAN MIX AND MATCH SHAPES TO EQUAL OR LESSER VALUE.
 NO CHARGE PRODUCT TO EQUAL OR LESSER VALUE.
 VALID FOR 3, 5, AND 25 PACKS OFFER EXCLUDES STANDARD CARBIDES, GREAT WHITE® GOLD AND TRIMMING AND FINISHING BURS.



STERILE R AVAILABLE IN ROLLS OF 25

Piranha® - Reduce sterilization costs & increase office productivity with a fresh diamond for each patient

Jazz® Polishers - Create faster, more beautiful polishes and eliminate up to 3 steps when finishing and polishing



All Lab Products

BUY 4 GET 1 AT NO CHARGE**

CHOOSE FROM: GREAT WHITE® LAB CUTTERS, JAZZ® LAB POLISHERS, EXPRESS LINE AND EXPRESS LINE-H DIAMOND DISCS, DURADI™ AND LAZER SHARP™ CAD/CAM BURS

Can mix and match shapes to equal or lesser value - No charge product to equal or lesser value

©2023 SS White® All Rights Reserved. SS White®, Express Line®, Great White® and Jazz® are registered trademarks of SS White Dental. DuraDi™ and Lazer Sharp™ are trademarks of SS White Dental. *Market research and test studies available upon request.**No charge goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. *SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim no charge goods: Mail a copy of your invoice by January 31, 2024 to: SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or email freegoods@sswhitedental.com. Please note that these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery. SS White® does its best to ensure that its photos/illustrations are as true to color as possible. However, due to inconsistencies of various monitors, lighting sources, digital photography, and manufacturing lot variations, SS White cannot guarantee that the color that appears on your screen/printed material accurately portrays the true color of the received product. Images are intended as a guide only.



Better Patient Outcomes
 Improved Efficiency
 Faster Practice Growth

1145 Towbin Avenue Lakewood, New Jersey 08701
 Visit us on the web at www.sswhitedental.com

Clean, Safe and Efficient. SS White® Oral Surgery Burs

SS White® Oral Surgery Burs are packaged sterile to save pre-operative time for the oral surgery practice. SS White® Oral Surgery Burs have been shown to save 3-5 minutes per procedure* by combining both the removal of bone and sectioning impacted third molars in one step, utilizing the popular 1703 bur which combines the sculpting ability of a round bur with a high performance 703 sectioning bur.

STERILE R ALL ORAL SURGERY BURS AVAILABLE IN PACKS OF 10



**BUY 10 PACKS
GET 2 PACKS
AT NO CHARGE****

**BUY 30 PACKS
GET 10 PACKS
AT NO CHARGE****

ORAL SURGERY BUR SHAPE SELECTION



Round		NON-STERILE Shank #1		STERILE Shank #2	STERILE Shank #3	STERILE Shank #4		STERILE Shank #5		STERILE Shank #6	
Size	Head Dia.	SSW#	Sterile Package*	SSW#	SSW#	SSW#	STRYKER#	SSW#	STRYKER#	SSW#	STRYKER#
1	1.0	14823									
2	1.2	14826									
4	1.4	14832									
6	1.8	14838	34838	30020	30035	30050	2296-100-418	30060	2296-101-418	30000	277-10-418
8	2.3	14844	34844	30021	30036	30051	2296-100-423	30031	2296-101-423	30001	277-10-423
10	2.7	14846	64846								



Cross Cut Fissure		Shank #1		Shank #2	Shank #3	Shank #4		Shank #5		Shank #6	
Size	Head Dia.	SSW#	Sterile Package*	SSW#	SSW#	SSW#	STRYKER#	SSW#	STRYKER#	SSW#	STRYKER#
557	1.0	14871									
558	1.2	14874									
560	1.6	14880				30052				30002	
1557	1.0	14898									
1558	1.2	14901									
700	1.2	14883									
701	1.2	14886	34916	30022		30053	2296-100-212	30062			
701L	1.2	14910									
702	1.6	14889	34889	30023	30037	30054	2296-100-216	30063	2296-101-216	30003	277-10-216
702L	1.6	14890		30024		30055	2296-100-316	30064	2296-101-316	30004	
1702	1.6	14916		30029		30048		30069			
1702L	1.6	14917		30030		30049		30070			
703	2.1	14892	34892	30025	30038	30056	2296-100-221	30065	2296-101-221	30005	277-10-221
703L	2.1	14893		30026	30039	30057	2296-100-321	30066	2296-101-321	30006	
1703	2.1	14914	34914	30027	30040	30058		30067		30007	
1703L	2.1	14915	34915	30028		30059	2296-100-421	30068	2296-101-421	30008	

Oral Surgery Bur Shank Length #1-6

Shank #1 - 44.5mm

Fits Hall® & Most Standard Hand Pieces for Medium Bur Guard



Shank #2 - 51mm

Fits Hall® & MicroAire® Hand Pieces with Medium Bur Guard



Shank #3 - 65mm

Fits Hall® & MicroAire® Hand Pieces with Long Bur Guard



Shank #4 - 44.5mm

Fits Stryker® 100K & Impaction Handpiece



Shank #5 - 59mm

Fits Stryker® 100K & Impaction Handpiece



Shank #6 - 44.5mm

Fits Stryker® 50K & Osteotome Drill J-Notch



**Can mix and match shapes to equal or lesser value - No charge product to equal or lesser value.

©2023 SS White® All Rights Reserved. SS White®, is a registered trademarks of SS White Dental. Hall® is a registered trademark of Zimmer, Inc. • Stryker® is a registered trademark of Stryker Corporation • MicroAire® is a registered trademark of MicroAire Surgical Instruments, Inc.*Market research and test studies available upon request.**No charge goods must be of equal or lesser value. All offers are limited to stock on hand and are not valid on prior purchases. Offers may not be combined with any other sale, promotion, discount, code, coupon and/or offer. Sales promotion choice must be indicated on invoice sent to SS White® on applicable order. SS White® reserves the right to discontinue a promotion at any time. †SS White's 100% guarantee is on all of our products. If you are not completely satisfied with your purchase, SS White® will issue a full refund of the purchase price through your distributor. To claim no charge goods: Mail a copy of your invoice by January 31, 2024 to: SS White Burs, Inc., Attn: Promotions Department, 1145 Towbin Ave., Lakewood, NJ 08701 or email freegoods@sswhitedental.com. Please note that these are end-user doctor promotions and must be redeemed by the qualifying practice. Please allow 4-6 weeks for free delivery. SS White® does its best to ensure that its photos/illustrations are as true to color as possible. However, due to inconsistencies of various monitors, lighting sources, digital photography, and manufacturing lot variations, SS White cannot guarantee that the color that appears on your screen/printed material accurately portrays the true color of the received product. Images are intended as a guide only.

Packaged Sterile

For your convenience and safety, SS White® Oral Surgery burs shank #1-6* come in a single sterile package.



100% Satisfaction Guaranteed†



Better Patient Outcomes
Improved Efficiency
Faster Practice Growth

1145 Towbin Avenue Lakewood, New Jersey 08701
Visit us on the web at www.sswhitedental.com